



Providing High-Value Care

Promote the COPDhub App

PROJECT OVERVIEW

STEP 1

Identify every patient on your COPD register.

STEP 2

Refer all patients to COPDhub app.

MORE DETAILS

STEP 1 - IDENTIFY YOUR COPD REGISTER

Identify every patient on your COPD register. The recommended codes for the diagnosis of COPD is:

- Read code: H3...00
- SNOMED code: 13645005

STEP 2 - PROMOTE THE APP

Using text messages, emails or printed leaflets (linked in the resources section) ensure the population of patients with COPD on your registers are aware of the COPDhub app as a mechanism for learning more about their COPD, checking their inhaler technique, monitoring their symptoms for signs of deterioration and a supplementary exercise programme that they can try at home.